



LEXCE Trend
















Exceptional processing unit with integrated drill

The drill unit uses a 5-axis mechanism, providing a high degree of accuracy for all your drilling jobs. The processing unit that runs the drill, also performs high quality safety beveling and grooving on any lenses.

- · 3D drilling optimally controlled by 5-axis
- · Multiple hole types covering extensive frames
- · Drilling angle can be set automatically or manually
- · Three types of drill bit (optional) for perfect fit
- · Precise grooving providing attractive edge surface regardless of lens shape







Intelligent blocker with integrated imager

Blocker unit is simple to operate while offering great performance. The integrated imager can capture optical tracings, along with drill hole data. The data can be easily edited on the multifunction color screen.

- \cdot Dual lens stage allows settings of all lens types
- $\cdot \ \text{Magnification of the display facilitates viewing of lens markings during blocking} \\$
- \cdot Highly accurate and precise blocking function
- \cdot Automatic hole and shape data acquisition by imager (optical tracer)
- · Screen enlargement facilitates hole data editing









NIDEK original design 3D frame tracer performs highly precise measurements. Additionally, two types of tracing methods are available for tracing demo lenses and patterns with the LEXCE Trend.

- · 3D frame tracing with full auto clamping (optional)
- · High curve frame measurement
- · Frame holder keeps frame in natural state during tracing
- · Reliable demo lens and pattern measurement by imager (optical tracer)
- · Demo lens and pattern tracing by Radius Measuring Unit in processing chamber







A 7-inch color LCD touch screen displays lens shape and layout in full scale. Processing conditions can be intuitively entered on the screen.

- · User preference of operation can be pre-set via software interface
 - Wizard mode; step-by-step operation, for beginners
- Professional mode; single screen operation, for experts
- · Uniquely designed, clearly visible icons
- · High resolution color LCD touch screen
- · Capacitive technology touch screen improves sensitivity







Proven high quality finishing

Thanks to avant-garde design and engineering innovations, the LEXCE Trend is technologically advanced, offering consistency and size accuracy while encompassing a faster cycle-time.

- \cdot Wider wheel capable of processing high Rx lenses
- · Full estimate soft processing mode controls axis shift
- · Water rinsing cycle keeps grinding chamber clean at all times
- · Customizable mini bevel is ideal for thin, metal eyewire frames
- · Lens edge polishing for flat and bevel lenses
- · Special wheel design for high base curve lens processing*
- · Multi bevel shapes to meet today's challenging eyewire frames*
- · Mini step bevel to grind an asymmetrical shelf-style rear bevel*

*Available for the type PLB-2R8










Enhanced user productivity

The LEXCE Trend is perfect for facilities with limited space.

Multiple functions with well-combined features, all in a compact footprint, improves productivity.

- \cdot Next job setup function
- \cdot Shape rotation adjustment function
- · Shape editing mode
- · Memory function for shape data management
- · Feature-rich compact design
- · Auto processing chamber door
- · Lit processing chamber for high visibility
- · Cooling water control knob



"A LEXCE" for everyone

The best option can be selected from several configurations depending on individual needs.







Complete with standalone configuration

- · Covers all necessary functions to make eyewear efficiently
- · Perfect for facilities where workspace limitations exist







Blends well with an external blocker with imager function

- · Parallel workflow for blocking or editing while processing
- · Selectable imager function depending on connected blocker







Blends well with an external blocker without imager function

- · Parallel workflow for blocking or editing while processing
- · Selectable imager function depending on connected blocker







Adoptable for high-volume processing

- · Capable of becoming part of a server software system
- · Provides flexibly to increase the number of lens edgers



For high curve and drilling, refer to the comparison chart.

| Minimum grinding size | Pliable cup (standard) W x H mm | | Mini cup (optional) W x H mm | | Nano cup (optional) W x H mm | |
|--------------------------|------------------------------------|--------------|---------------------------------|--------------|---------------------------------|--------------|
| | PLB-2R8 | PLB-2R | PLB-2R8 | PLB-2R | PLB-2R8 | PLB-2R |
| Flat edging | ø32.0 x 19.0 | ← | ø22.0 x 17.4 | ← | ø20.0 x 15.5 | \leftarrow |
| Bevel edging | ø33.0 x 20.6 | \leftarrow | ø23.0 x 18.4 | \leftarrow | ø21.0 x 16.5 | \leftarrow |
| Safety beveling (flat) | ø35.0 x 22.0 | \leftarrow | ø25.0 x 20.3 | \leftarrow | ø23.0 x 18.5 | \leftarrow |
| Safety beveling (bevel) | ø36.6 x 23.6 | \leftarrow | ø26.6 x 21.9 | \leftarrow | ø24.6 x 20.1 | \leftarrow |
| High base curve beveling | ø37.8 x 24.8 | | ø27.8 x 23.2 | | ø25.8 x 21.3 | |
| Grooving | ø32.0 x 19.0 | \leftarrow | ø22.0 x 17.4 | \leftarrow | ø20.0 x 15.5 | \leftarrow |

▶ Mate 1 with ICE mini+ style



▶ Mate 2 with Ice 900 style



System configurations

▶ **Pro** Combination style – Mini lab



▶ **Pro** Combination style – Extended lab



| Type comparis | son | | Frame | Shape imager | Blocker | Grooving | High curve | Drilling |
|-----------------------|----------|-----------|-------|--------------|---------|----------|------------|-----------------|
| PLB-2R8 Core PLB-2R | Drill | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | F LB-2NO | Non drill | ✓ | ✓ | ✓ | ✓ | ✓ | |
| | DLD 2D | Drill | ✓ | ✓ | ✓ | ✓ | | ✓ |
| | PLD-ZN | Non drill | ✓ | ✓ | ✓ | ✓ | | |
| PLB-2R8 | DI D 200 | Drill | ✓ | ✓ | | ✓ | ✓ | ✓ |
| | FLB-ZRO | Non drill | ✓ | ✓ | | ✓ | ✓ | |
| Mate 1 | | Drill | ✓ | ✓ | | ✓ | | ✓ |
| PLB-2 | PLB-2R | Non drill | ✓ | ✓ | | ✓ | | |
| PLB-2R8 Mate 2 PLB-2R | DI D 200 | Drill | ✓ | | | ✓ | ✓ | ✓ |
| | PLD-2NO | Non drill | ✓ | | | ✓ | ✓ | |
| | DLD 2D | Drill | ✓ | | | ✓ | | ✓ |
| | PLB-2R | Non drill | ✓ | | | ✓ | | |
| PLB-2F | DI D 200 | Drill | | | | ✓ | ✓ | ✓ |
| | PLD-ZNO | Non drill | | | | ✓ | ✓ | |
| PLB-2R | DLP 2P | Drill | | | | ✓ | | ✓ |
| | PLB-ZN | Non drill | | | | ✓ | | |

LEXCE Trend Specifications

| Model Grinding system | LEXCE Trend8 Patternless | LEXCE Trend ← | | |
|----------------------------------|--|---|--|--|
| Mode | Beveling (automatic, guided, safety beveling, polishing, high base curve), Flat edging (polishing, safety beveling, grooving), Drilling, Mini beveling (0.4 to 0.7 mm) (0.1 mm increments), Mini step processing (0.0 to 3.8 mm) (0.1 mm increments), Custom beveling, Frame changing, Soft processing | Beveling (automatic, guided, safety beveling, polishing), Flat edging (polishing, safety beveling, grooving), Drilling, Mini beveling (0.4 to 0.7 mm) (0.1 mm increments) Frame changing, Soft processing | | |
| letting range | custom bevening, frame changing, sort processing | | | |
| FPD | 30.00 to 99.50 mm (0.01 mm increments) | | | |
| PD | 30.00 to 99.50 mm (0.01 mm increments) | | | |
| 1/2PD | 15.00 to 49.75 mm (0.01 mm increments) | ← | | |
| Optical center height | 0 to ±15.0 mm (0.1 mm increments) | | | |
| Size adjustment | 0 to ±9.95 mm (0.01 mm increments) | | | |
| Bevel position | 0 to ±9.95 mm (0.01 mm increments) | | | |
| Ninimum grinding size | | | | |
| Flat edging | ø32.0 x 19.0 mm / with nano cup (optional) ø20.0 x 15.5 mm | ← | | |
| Bevel edging | ø33.0 x 20.6 mm / with nano cup (optional) ø21.0 x 16.5 mm | ← | | |
| Safety beveling (flat) | ø35.0 x 22.0 mm / with nano cup (optional) ø23.0 x 18.5 mm | ← | | |
| Safety beveling (bevel) | ø36.6 x 23.6 mm / with nano cup (optional) ø24.6 x 20.1 mm | ← | | |
| High base curve beveling | ø37.8 x 24.8 mm / with nano cup (optional) ø25.8 x 21.3 mm | None | | |
| Grooving | ø32.0 x 19.0 mm / with nano cup (optional) ø20.0 x 15.5 mm | ← | | |
| Drilling*1 | | | | |
| Hole diameter | ø0.80 to 10.00 mm (0.01 mm increments) | | | |
| Hole depth | 6.0 mm or less | | | |
| Range for hole milling | ø34.0 to 68.5 mm from lens rotation axis | | | |
| Direction for hole milling | Automatic/Manual tilting 2.5 to 18° | ← | | |
| Slotted hole width | ø0.80 to 10.00 mm (0.01 mm increments) | | | |
| Slotted hole depth | 6 mm or less | | | |
| Slotted hole length | 20 mm or less | | | |
| Blocking unit*2 | | | | |
| Method | Manual blocking | ← | | |
| Blocking position accuracy | ±0.5 mm | | | |
| Axis angle accuracy | ±1.0° | | | |
| hape imager function*3 | | | | |
| Measurement range | 65.0 x 50.0 mm (±1.5 mm) | | | |
| Hole position | 0.01 mm increments | ← | | |
| Hole diameter | ø0.80 to 10.00 mm (0.01 mm increments) | | | |
| Demo lens / pattern tracing unit | | | | |
| Method | Shape measurement using feeler unit | | | |
| Measuring points | 1,000 points | ← | | |
| Measurement range | ø22.0 to 76.0 mm (17.4 to 66.0 mm vertically) | | | |
| rame tracer (optional) | | | | |
| Method | Automatic 3D binocular tracing | | | |
| Measuring points | 1,000 points | | | |
| Measurement range | Shape width : 23.0 to 70.0 mm | | | |
| | Shape height : 18.4 to 66.0 mm | ← | | |
| | Frame horizontal width: 113 to 150 mm | , · | | |
| FPD measurement | Available | | | |
| Frame clamping | One-touch automatic clamping | | | |
| Setting of stylus | Switchable between automatic and semiautomatic | | | |
| Measurement accuracy | Frame tracing ±0.1 mm | | | |
| Wheel configuration | Type PLB-2R8 | Type PLB-2R | | |
| Vater supply system | Pump circulation or direct connection to tap water | ← | | |
| nterface | RS-232C - 1 port | | | |
| | LAN -1 port | ← | | |
| | USB - 1 port | , | | |
| Power supply | 100 to 120 V / 240 V AC, 50/60 Hz | ← | | |
| Power consumption | 1.3 kVA | - | | |
| Dimensions/Mass | 545 (W) x 530 (D) x 460 (H) mm / 38.5 kg (Core), 37.8 kg (Mate 1) | | | |
| | 545 (W) x 434 (D) x 460 (H) mm / 37.2 kg (Mate 2) | | | |
| | 545 (W) x 434 (D) x 344 (H) mm / 34.6 kg (Pro) | ← | | |
| | 21.5 (W) x 20.9 (D) x 18.1 (H)" / 84.9 lbs. (Core), 83.3 lbs. (Mate 1) | | | |
| | 21.5 (W) x 17.1 (D) x 18.1 (H)" / 82.0 lbs. (Mate 2) | | | |
| tandard accessories | 21.5 (W) x 17.1 (D) x 13.5 (H)" / 76.3 lbs. (Pro) Drill bit (10 units)*1, Hexagonal screwdriver (2.5 mm), Hexagonal wrench (2.0 mm, 3.0 mm, | | | |
| tarruaru accessories | and 4.0 mm), Dressing stick for glass roughing wheel, Dressing stick for finishing wheel, | | | |
| | Compound kit for polishing wheel, Pliable cup, Pliable cup for high base curve lenses, | ← | | |
| | | <u></u> | | |
| | | | | |
| | Double-coated adhesive pad, Pliable cup remover, Adapter set, Pattern holder, Stage for | | | |
| Ontional accessories | Double-coated adhesive pad, Pliable cup remover, Adapter set, Pattern holder, Stage for small diameter lens*3, Calibration jig, Flat lens, Ferrite core, Accessory case, Power cord | | | |
| Optional accessories | Double-coated adhesive pad, Pliable cup remover, Adapter set, Pattern holder, Stage for | ← | | |

^{*1} Available for the drill-equipped model

Specifications and design are subject to change without notice.



HEAD OFFICE (International Div.) 34-14 Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038, JAPAN TEL: +81-533-67-8895 URL: www.nidek.com [Manufacturer]

TOKYO OFFICE (International Div.) 3F Sumitomo Fudosan Hongo Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN TEL: +81-3-5844-2641

URL: www.nidek.com

NIDEK INC. 2040 Corporate Court, San Jose, CA 95131, U.S.A. TEL: +1-408-468-6400 +1-800-223-9044 (US Only) URL: usa.nidek.com

NIDEK S.A. Europarc, 13 rue Auguste Perret, 94042 Créteil, FRANCE TEL: +33-1-49 80 97 97 URL: www.nidek.fr

Via dell'Artigianato, 6/A, 35020 Albignasego (Padova), ITALY

TEL: +39 049 8629200/8626399 URL: www.nidektechnologies.it

NIDEK TECHNOLOGIES S.R.L. NIDEK (SHANGHAI) CO., LTD. NIDEK SINGAPORE PTE. LTD. Rm3205,Shanghai Multi Media Park, No.1027 Chang Ning Rd, Chang Ning District, Shanghai, CHINA 200050 TEL: +86 021-5212-7942 URL: www.nidek-china.cn

51 Changi Business Park Central 2, #06-14, The Signature 486066, SINGAPORE TEL: +65 6588 0389 URL: www.nidek.sg

^{*2} Available for the Core

^{*3} Available for the Core and Mate 1