



Specular Microscope
CEM-530



THE ART OF EYE CARE





CEM-530

SPECULAR MICROSCOPE

Paracentral Specular Microscopy

In addition to conventional central and peripheral specular microscopy, the CEM-530 includes a unique function to capture paracentral images. The paracentral images are captured at eight points, 5° visual angle within a 0.25 mm x 0.55 mm field and enable enhanced assessment surrounding the central image.




Simulated image of 15 fixation lights*
Central 1 point
Paracentral 8 points (5° visual angle)
Peripheral 6 points (27° visual angle)
*Only one selected fixation light is on.

Paracentral mode provides a total image of endothelial cells.

The paracentral mode allows detailed evaluation of cell shape, which is important for preoperative assessment. For example, assessment of corneal guttata using a central image only is often clinically ineffective due to the limited number of countable cells.

Supervisor: Prof. Yuichi Ohashi
Department of Ophthalmology,
Ehime University School of Medicine



Paracentral Image

Faster Measurements and Two-second Auto Analysis

The CEM-530 with new advanced software and enhanced image capture system allows rapid image acquisition. A new algorithm for the automated analysis software performs complete analysis in two seconds.



Comprehensive Analysis

The analysis results with graphic and color-coded cell images helps the clinician to rapidly and effectively evaluate the endothelial cell layer.




Detail analysis

Advanced Manual Analysis Functions

Two new methods, center point and corner point, have been added to the manual analysis function. These additions provide the clinician with three manual analysis functions.


Center point

Select the approximate center of a cell. The cells are detected based on the surrounding points. This method is effective for areas where groups of cells are clumped together.




Corner point

Trace the outlines of the cells to be analyzed by selecting the corners of each cell. This method is suitable for detailed identification of the size and dimension of isolated cells.



Pattern select

Select a hexagonal reference pattern that is similar to the cell size and drag it onto the cell to be analyzed. This method is effective for rough identification of the size and dimension of the cells.




Combination of Auto and Manual Analyses

All three manual analysis methods can be performed on the same image and also on auto-analyzed images. The versatility of combining automated and manual analyses allows analysis of the range of pathology in a comprehensive practice.

Corner point

Auto analysis


Center point



Easy Operation

3-D Auto Tracking, Auto Shot, and Tiltable Touch Screen

The 3-D auto tracking, auto shot, and tiltable touch screen provide ease of use, allowing faster and more accurate measurement.



Instant Printout with Built-in Printer


The built-in printer provides an instant printout of the analyzed data and images of the endothelial cells.



CEM Viewer for NAVIS-EX

CEM Viewer is software for viewing and working with CEM-530 data within NAVIS-EX, image filing software of NIDEK. This functionality enhances the capability of the CEM-530 with additional features and increases clinical efficiency.

- Data Management and Endothelial Cell Count
- Progression Follow-up and Comparison
- Paracentral Display with Peripheral



CEM-530 Specifications

Endothelial image capture	
Capture field	0.25 (W) x 0.55 (H) mm
Capture position	Central 1 point Paracentral 8 points (5° visual angle, 45° spacing) Peripheral 6 points (27° visual angle, 60° spacing)
Pachymetry	
Measurement range	300 to 1,000 µm
Accuracy	±10 µm
Auto tracking	X-Y-Z directions
Auto shot	Available
Display	Tilttable 8.4-inch color LCD touch screen
Printer	Built-in thermal line printer External video printer (optional)
Interface	LAN, USB, Video output (BNC connector for video printer)
Power supply	AC 100 to 240 V 50 / 60 Hz
Power consumption	100 VA
Dimensions / Mass	291 (W) x 495 (D) x 457 (H) mm / 20 kg 11.5 (W) x 19.5 (D) x 18.0 (H) " / 44 lbs.



Product / Model name: Specular Microscope CEM-530

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.



HEAD OFFICE
(International Div.)
34-14 Maehama, Hiroishi
Gamagori, Aichi 443-0038,
JAPAN
TEL: +81-533-67-8895
URL: <http://www.nidek.com>
[Manufacturer]

TOKYO OFFICE
(International Div.)
3F Sumitomo Fudosan Hongo
Bldg., 3-22-5 Hongo, Bunkyo-ku,
Tokyo 113-0033, JAPAN
TEL: +81-3-5844-2641
URL: <http://www.nidek.com>

NIDEK INC.
47651 Westinghouse Drive,
Fremont, CA 94539, U.S.A.
TEL: +1-510-226-5700
+1-800-223-9044
(US only)
URL: <http://usa.nidek.com>

NIDEK S.A.
Europarc,
13 rue Auguste Perret,
94042 Créteil, FRANCE
TEL: +33-1-49 80 97 97
URL: <http://www.nidek.fr>

NIDEK TECHNOLOGIES S.R.L.
Via dell'Artigianato,
6/A, 35020 Albignasego (Padova),
ITALY
TEL: +39 049 8629200/8626399
URL: <http://www.nidektechnologies.it>

NIDEK (SHANGHAI) CO., LTD.
#915, China Venturetech Plaza,
819 Nanjing West Rd, Jing An
District, Shanghai 200041,
CHINA
TEL: +86 021-5212-7942
URL: <http://www.nidek-china.cn>

NIDEK SINGAPORE PTE. LTD.
51 Changi Business Park
Central 2, #06-14, The
Signature 486066,
SINGAPORE
TEL: +65 6588 0389

