


# SC-1600 (Pola) / 1600



THE ART OF EYE CARE

# Reliable Features

----

for easy operation and efficient workflow

# Sufficient Functions

for more refined and smoother examinations

# Flexible Usability

for comfortable test in any circumstances



## SC-1600@ / 1600 Chart Types

	0.2.0       0.04 to 2.0         0.1 to 1.25         0.1 to 1.2            1.25       0.1 to 1.0         0.2.0       0.32 to 2.0	Type M*1         500 to 10         200 to 10         100 to 10         —         200 to 20         64 to 10         —         @	Type UK*1         150 to 4         12 to 4            150 to 4         38 to 6         20 to 3         Image: Comparison of the second sec	Type G+*1 0.04 to 1.6 0.04 to 1.6 0.04 to 1.6 0.04 to 1.6 0.1 to 1.6 0.32 to 2.0	Type U*3         400 to 10         200 to 10         100 to 15            100 to 20            ioi:         @	Type F *3         0.05 to 1.2         0.2 to 1.2         0.1 to 1.4            ***         •         •         •         •         •	Type ISO*3
0 1.5       —         0 1.5       0.04 to 1         0 1.0       0.1 to 1         0 2.0       0.32 to 2	0.1 to 1.25 0.1.6 0.1 to 1.2 	200 to 10 100 to 10 	12 to 4 — 150 to 4 38 to 6 20 to 3 … … … … … … … … … … … … …	0.04 to 1.6 0.04 to 1.6 0.04 to 1.6 0.1 to 1.6 0.32 to 2.0 	200 to 10 100 to 15 — 100 to 20 — •••••••••••••••••••••••••••••••••••	0.2 to 1.2 0.1 to 1.2 0.1 to 1.2 0.1 to 1.0 — — — — — — —	0.1 to 1.25 0.1 to 1.6 0.16 to 1.0 — — — —
1.5       0.04 to 1         .	0.1.6       0.1 to 1.2             1.25       0.1 to 1.0         0.2.0       0.32 to 2.0         4          •       •	100 to 10 		0.04 to 1.6 0.04 to 1.6 0.1 to 1.6 0.32 to 2.0 	100 to 15 — 100 to 20 —	0.1 to 1.2 0.1 to 1.2 0.1 to 1.0 — — — — — —	0.1 to 1.25 0.1 to 1.6 0.16 to 1.0 — — — —
			38 to 6 20 to 3 	0.04 to 1.6 0.1 to 1.6 0.32 to 2.0 	 100 to 20  @	0.1 to 1.2 0.1 to 1.0 — — — — —	0.1 to 1.6 0.16 to 1.0 — — —
0.32 to 2	0.32 to 2.0       0.32 to 2.0       • <td< td=""><td>64 to 10</td><td>38 to 6 20 to 3 </td><td>0.1 to 1.6 0.32 to 2.0</td><td></td><td>0.1 to 1.0</td><td>0.16 to 1.0</td></td<>	64 to 10	38 to 6 20 to 3 	0.1 to 1.6 0.32 to 2.0		0.1 to 1.0	0.16 to 1.0
0.32 to 2	0.32 to 2.0       0.32 to 2.0       • <td< td=""><td>64 to 10</td><td>20 to 3</td><td>0.32 to 2.0</td><td></td><td>- **</td><td>- *</td></td<>	64 to 10	20 to 3	0.32 to 2.0		- **	- *
			Image: Constraint of the second sec	**** * *	@	** •	۲
			* ***	* ****		•	۲
		_					
•••••	·····				—	-	
		*****				1	
				<u>+</u>			<u></u>
					—	10.0X	
+	_	_	_	÷ +	$\div$ +	+ +	_
·			_				
	_		_	_		_	_
	_	*****	_	_		_	_
. –			<b>–</b>	•	•	•	•
		_	_	_	_	•	_
		_	_	(+)	_	+	Ŧ
		• 11 • 11 • 11 • 11			• 11 • 11 •		
		+:+	+++	•••	•••	+;+	+:-
Mallet	et phoria) )	_	H K K	Pointer Double pointer II Stereo Stereo triangle balance	_	_	_
	Mallı (horizontal ()	Mallet (horizontal phoria)	Image: Constraint of the second se	Image: Second	Image: Constraint of the constr	Image: set of the set of	Image: Second

\*1 only with the RT-5100 / 3100  $\,$  \*2 only for the SC-1600Pola  $\,$  \*3 only with the RT-2100 / 900  $\,$ 

## Wide Range of Image Display

The high definition 17-inch LCD screen allows VA 0.04 to 2.0 images to be displayed.

The wide bezel secures enough space around

the chart to reduce peripheral fusion. The microcomputer-based chart operates quickly just

after start up.



#### Programmed Test

For easy refraction, a maximum of three test programs can be set to the remote control according to the user's preference.

#### Auto Light Off Mode

For power saving, the screen backlight can be set to go out automatically after the unit is idle for 15 minutes.

#### Contrast Test

To support patients who have undergone cataract or refractive surgery, the SC-1600Pola / 1600 has the ability to measure contrast sensitivity at three different levels below the normal threshold.



#### **Center Display Function**

Single letter and horizontal / vertical line masking image are always centered on the screen. Characters displayed in a horizontal line can be randomized to prevent letter memorization.



3 E M E W
l line






#### Black and White Inversion Function

The VA chart can be selected from black on white to white on black.

#### Selectable Refraction Distance

The refraction distance can be pre-selected from 2.5 to 6 m in 1 cm or 1 inch increments to fit any size of exam room. The chart image size can be reduced or enlarged corresponding to the refraction distance to keep the same visual angle. The same refraction distance is also adjustable to accommodate mirrored exam rooms.



## System Integration

The NIDEK system chart can be synchronized with the NIDEK Intelligent Refractor by optional connection cable or infrared communication unit for more refined

W ir p Int con

and smoother examinations. When synchronized, chart image adjustment is possible using the Intelligent Refractor control box.

#### Auto Test Mode

By simply pressing Yes / No button on the remote control in accordance with the patient response, the test proceeds

automatically with no need to select the next chart.



#### Indication of a Slide Image

When not being used for eye examination, Images stored in an SD memory card can be displayed. It is also possible to advertise the PR of a shop and the article by the slide function. There are a lot of 3-D images for display.




#### Night Vision Test

The NIDEK SC-1600Pola / 1600 provides "Night Vision Test" to examine night-time visual acuity under low light condition.



#### Polarized Charts for Binocular Test\*

The SC-1600Pola adopts the state-of-the-art technology, special LCD to polarize charts, which provides binocular assessment in a conventional way with stylish design.





\*only for the SC-1600Pola

#### Stand Options

The optional floor stand and table stand expand the utility.



#### Operation with Remote Control

The SC-1600 series can be operated with a newly designed optional remote control. Its ergonomic design makes it easy to handle, and the large LCD has room to display plenty of information. A simple but well-organized button layout provides intuitive operability.



#### SC-1600 @ / 1600 Specifications

Model	SC-1600@	SC-1600	
LCD			
Туре	17-inch SXGA liquid crystal display	←	
Resolution	1280 x 1024 pixels		
Brightness	200 cd / m <sup>2</sup> (night mode 10 cd / m <sup>2</sup> )		
Chart type	T, PhM, P, M, UK, G+, U, F, ISO	T, PhM, P, M, UK	
Distance of installation	2.5 to 6 m, 1 cm / 1" increments	←	
Binocular vision test	Polarized chart	Red-green, Prism	
Chart changeover devices	Remote control (optional)	←	
	Control box for the NIDEK Intelligent Refractor RT-5100, RT-3100		
Memory card	SD memory card / Built-in card slot	←	
Power supply	AC 100 to 240 V	←	
	50 / 60 Hz	N N N N N N N N N N N N N N N N N N N	
Power consumption	100 VA	←	
Dimensions / Mass	520 (W) x 47 (D) x 455 (H) mm / 5 kg	←	
	20.5 (W) x 1.9 (D) x 17.9 (H) " / 11.0 lbs.	```	
Standard accessories	Polarizing glasses, Tilt bracket, Screws, Hanging plate, Wood screws,	Hanging plate, Wood screws, Power cord, Positioning paper	
	Power cord, Positioning paper		
Optional accessories	Red-green glasses, Reflection mirror, Floor stand,	Red-green glasses, Prism glasses, Reflection mirror, Floor stand,	
	Table stand, Communication cable for the refractor,	Table stand, Communication cable for the refractor,	
	Infrared communication unit for the refractor, Remote control	Infrared communication unit for the refractor, Remote control	

Product / Model name: SYSTEM CHART SC-1600 Brochure and listed features of the device are intended for non-US practitioners. Specifications may vary depending on circumstances in each country. Specifications and design are subject to change without notice.



#### HEAD OFFICE (International Div.) 34-14 Maehama, Hiroishi Gamagori, Aichi 443-0038, JAPAN TEL: +81-533-67-8895 URL: http://www.nidek.com [Manufacturer]

#### TOKYO OFFICE (International Div.) 3F Sumitomo Fudosan Hongo Bldg, 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN TEL: +81-3-5844-2641

URL: http://www.nidek.com

#### NIDEK INC. 47651 Westinghouse Drive, Fremont, CA 94539, U.S.A. TEL:+1-510-226-5700 +1-800-223-9044 (US only) URL: http://usa.nidek.com

 44539, U.S.A.
 13 rue Auguste Perret,

 26-5700
 94042 Créteil, FRANCE

 23-9044
 TEL: +33-1-49 80 97 97

 URL: http://www.nidek.fr

NIDEK S.A.

Europarc,

NIDEK TECHNOLOGIES S.R.L.NIDEK (SHANGHAI) CO., LTD.NIDEK SINGAPORE PTE. LTD.Via dell'Artigianato,#915, China Venturetech Plaza,51 Changi Business Park6/A, 35020 Albignasego (Padova),819 Nanjing West Rd, Jing AnCentral 2, #06-14, The

94042 Créteil, FRANCE ITALY TEL: +33-1-49 80 97 97 TEL: +39 049 8629200/8626399 URL: http://www.nidek.fr URL: http://www.nidektchnologies.it

NIDEK (SHANGHAI) CO., LTD. #915, China Venturetech Plaza, 819 Nanjing West Rd, Jing An District, Shanghai 200041, CHINA TEL: +86 021-5212-7942 URL: http://www.nidek-china.cn NIDEK SINGAPORE PTE. LTD 51 Changi Business Park Central 2, #06-14, The Signature 486066, SINGAPORE TEL:+65 6588 0389

