

# US-4000 / 500



THE ART OF EYE CARE



## Three in One\*

The US-4000 comprehensive ultrasound device includes B-scan, Biometry and Pachymetry. Combining all three functions in one compact and innovative device brings a new image to your practice.

\*Only for the US-4000 The US-500 includes Biometry and Pachymetry.

•ONLY 8.5 kg COMPACT

- TOUCH PANEL DISPLAY NO PC REQUIRED
- TILTABLE COLOR LCD



## B-scan\*

## A Crisp and Clear Image

The results are high quality images, which are essential for accurate analysis. The key to achieving high quality images is through scanning 400 lines over 60° which is displayed on a 1024 x 768 XGA monitor. High contrast of resonated signals and the high quality LCD monitor provide the ideal images for the most accurate diagnosis.







Cross Vector Mode

### Internal Printer

B-scan image is easily printed with the internal thermal printer and may be attached to the patient's medical record sheet.



## Biometry

### Rapid, Accurate and Easy to Use

Using new algorithms, axial length measurements and IOL power calculations are performed twice as rapidly as the conventional model. Measurement accuracy for dense cataracts and existing opacities is accomplished by pressing the "Dense Cat" switch and "Gate switch".

Axial: 24.	45 mm ACC		75 mm			21 mm	Eye Type
No.10 Eye Type:Phakic Mode :Auto	Th. level :Normal Gain : 70		Asial V. 1550m/s Lens V. 1641m/s ACD V. 1532m/s			Phakic	
	Dense Cat:OFF	No.	Axial	ACD	Lens	Vit	Mode
		1	24.45	3.75	4. 21	16.49	Auto
		2	24.45	3.75	4.21	16.49	Auto
		3	24. 45	3.75	4, 21	16.49	Dense Cat
		4	24, 45	3.75	4, 21	16.49	Dense Cat
		5	24.45	3.75	4, 21	16.49	OFF
·		6	24.45	3.75	4, 21	16.49	
		7	24. 45	3.75	4, 21	16.49	Threshold
			24.45	3.75	4.21	16.49	Normal
		9	24.45	3.75	4, 21	16.49	HOTHER
		* 410	24, 45	3,75	4, 21	16.49	- Knob -
		* Avg	24.45	0.00	4, 21	0.00	Knob 1:
$\Delta \Delta \Delta$		05/L	0.00	0.00	0.00	0.00	Gate L/F
Retina Gate		Sele	et	Del	ete	Clear	Knob 2: Gain: 70

**Biometry Image** 



**Customized IOL Calculation** 

## Pachymetry\*1

The US-4000 / 500 pachymetry mode offers precise measurement of the corneal thickness within a  $\pm 5 \ \mu m$  error.



45° fixed Probe\*2



45° detachable Probe



Straight Probe

\*1 The model without pachymetry is also available.

\*<sup>2</sup> Provided as standard accessories depending on the configuration.

## Easy Data Transfer

Data is easily stored through USB memory and LAN. Utilizing the special interface cable, keratometry values measured by any NIDEK ARK are automatically imported into the US-4000 / 500 for easy IOL power calculation.

#### US-4000 / 500 Specifications

Model	US-4000	US-500
B-scan		
Probe	10 MHz transducer, 10 frames / second	
Scan angle	60°	
Scan depth	Normal (35 mm / 1550 m/s), Long (50 mm / 1550 m/s)	
Sector line density	400 Lines	N
Zoom	x2.5, x5.0	Not available
Moving image record	Approx. 20 seconds	
Scale	Color, Gray scale 256 levels	
Gain / TGC	0 to 90 dB variable / 0 to -20 dB variable	
Gain curves	Log, Liner, S-curves	
Biometry		
Probe	10 MHz solid probe	
Internal fixation	LED (red)	
Measurement value	Axial length, Anterior chamber depth, Lens thickness, Vitreous	
	body length	←
Accuracy	0.1 mm	
Range	12 to 40 mm	
Minimum calculation step	0.01 D	
Built-in IOL formula	BINKHORST, HOLLADAY, SRK, SRK II, SRK/T, HOFFR Q	
Pachymetry*1		
Probe	10 MHz solid probe	
Accuracy	5 μm	
Range	200 to 1300 μm	←
Minimum indicated unit	1 μm	
Measured part	Corneal thickness up to 25 points can be memorized.	
IOP correction	Available	
Display	8.4-inch TFT color LCD (XGA: 1024 x 768)	→
Printer	Thermal type line printer (Easy loading and auto cutter)	
Interface	USB memory (1.1), LAN, RS-232C for KM communication,	USB memory (1.1), LAN, RS-232C for KM communication
	Video out (NTSC)	
Power supply	AC 100 to 120 V 10%, 230 V 10%, 50 / 60 Hz	→
Power consumption	70 VA	→
Dimensions / Mass	300 (W) x 285 (D) x 330 (H) mm / 8.5 kg	←
Standard accessories	Stylus, B-scan probe, A-scan probe (14610-E310), Foot switch, Test	Stylus, A-scan probe (14610-E310), Foot switch, Test piece (for
	piece (for Biometry measurement), Printer paper, Power cord,	Biometry measurement), Printer paper, Power cord, Dust cover,
	Ultrasonic gel, Dust cover, Spare fuse, Probe rest	Spare fuse, Probe rest
Optional accessories	Video printer, Video printer paper, Pachymetry probe (45° fixed	Pachymetry probe (45° fixed type)*, Pachymetry probe (45°
	type)*, Pachymetry probe (45° detachable type), Pachymetry probe	detachable type), Pachymetry probe (straight type), Test piece for
	(straight type), Test piece for 45° pachymetry probes*, Test piece for	45° pachymetry probes*, Test piece for straight pachymetry probe
	straight pachymetry probe, Probe stand, Barcode reader, Magnetic	Probe stand, Barcode reader, Magnetic card reader, IOL calculation
	card reader, IOL calculation formula (Hagis), RS-232C cable	formula (Hagis), RS-232C cable
	*Provided as standard accessories depending on the configuration.	*Provided as standard accessories depending on the configuration

\*1 The model without pachymetry is also available.

Caution : U.S. Federal Law restricts this device to sale, distribution and use by or on the order of a physician or other licensed eye care practitioner. Specifications may vary depending on circumstances in each country. Specifications and design are subject to change without notice.



#### HEAD OFFICE

34-14 Maehama, Hiroishi Gamagori, Aichi 443-0038, Japan Telephone :+81-533-67-6611 Facsimile :+81-533-67-6610 URL : http://www.nidek.co.jp

#### [Manufacturer]



#### TOKYO OFFICE (International Div.) 3F Sumitomo Fudosan Hongo Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan

3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan Telephone: +81-3-5844-2641 Facsimile: :+81-3-5844-2642 URL: http://www.nidek.com

#### NIDEK INC.

47651 Westinghouse Drive Fremont, CA 94539, U.S.A. Telephone : +1-510-226-5700 : +1-800-223-9044 (US only) Facsimile : +1-510-226-5750 URL : http://usa.nidek.com

#### NIDEK S.A.

Europarc 13, rue Auguste Perret 94042 Créteil, France Telephone :+33-1-49 80 97 97 Facsimile :+33-1-49 80 32 08 URL : http://www.nidek.fr

#### NIDEK TECHNOLOGIES Srl

Via dell'Artigianato, 6 / A 35020 Albignasego (Padova), Italy Telephone: +39 049 862200 / 8626399 Facsimile: +39 049 8626824 URL: http://www.nidektechnologies.it